


# Left main bifurcation treated with Biomatrix Alpha.

Fco. Javier Lacunza Ruiz.  
Interventional Cardiology Department.  
Hospital Universitario Virgen de la Arrixaca.  
Murcia. Spain.

# Left main bifurcation treated with Biomatrix Alpha.

- **Clinical data:**
  - 74 y.o. lady.
  - CVRF: Hypertension, hyperlipidemia.
  - Referred to ER for a 2 hours chest pain.
  - EKG compatible with posterolateral STEMI.






# Left main bifurcation treated with Biomatrix Alpha.

- Clinical data (II)
  - Right Angiogram:
    - No significant stenosis.
  - Left Angiogram (right radial artery):
    - Acute occlusion of proximal circumflex (CX).
    - Severe stenosis in distal left main, affecting the ostium of LAD & CX.


# Left main bifurcation treated with Biomatrix Alpha.




➤ Diffuse disease in RCA.



# Left main bifurcation treated with Biomatrix Alpha.




- Critical calcified lesion in ostial LAD.
- Diffuse disease in prox-mid LAD.





# Left main bifurcation treated with Biomatrix Alpha.



- Severe calcified lesion in distal LM.
- Acute occlusion of proximal CX.
- Severe lesion in ostial CX .




# Left main bifurcation treated with Biomatrix Alpha.

- Initial strategy:
  1. Restore flow in proximal CX.
  2. LM treatment
    1. Double stent strategy due to the severe stenosis in both ostium (LAD & CX).
- Change to 7French femoral approach




# Left main bifurcation treated with Biomatrix Alpha.

- PCI in the CX lesion.
  - EBU 3.5 7F. 2 Sion guidewire to LAD & CX.



➤ CX opens with the pass of the wireguide.



➤ High burden of thrombus inside the artery. Reopro infusion started.






# Left main bifurcation treated with Biomatrix Alpha.


- Predilatation of CX lesion and ostial CX with 2.0 & 2.5 mm balloon.



# Left main bifurcation treated with Biomatrix Alpha.




- TIMI 3 flow in distal CX.
- Good expansion of the ostial balloon.
- Next step: Stent to proximal CX.





# Left main bifurcation treated with Biomatrix Alpha.





- Result after a 3.0 x 29 mm *Biomatrix Alpha* stent implanted in prox CX.
- Good expansion of the stent. Distal TIMI 3 flow.
- Next step: Stent to LM with “Crushing stent” technique.





# Left main bifurcation treated with Biomatrix Alpha.

➤ “Crushing stent” technique:

➤ Predilatation of ostial LAD with 3.0 mm NC balloon.  
Good expansion.



# Left main bifurcation treated with Biomatrix Alpha.



➤ “Crushing stent” technique:

➤ Stent position:


**CX stent** (*Biomatrix Alpha* 3.0 x 12mm) overlapped with the previous stent and minimally prolapsing into LM.

**LAD stent** (*Biomatrix Alpha* 3.0 x 19 mm) from proximal LM to proximal LAD.

➤ Initial inflation of the CX stent.



# Left main bifurcation treated with Biomatrix Alpha.




➤ “Crushing stent” technique:

➤ Initial result after the implant of both stents:

TIMI 3 flow in both arteries  
with minimal residual  
stenosis.



# Left main bifurcation treated with Biomatrix Alpha.




➤ “Crushing stent” technique:

➤ Optimization of the results:

Recrossing towards CX and dilatation of the struts of the stent with 1.5 and 2.0 mm NC balloon.

# Left main bifurcation treated with Biomatrix Alpha.




➤ “Crushing stent” technique:

➤ Optimization of the results:

“Kissing balloon” with two 3.0 NC balloon.



# Left main bifurcation treated with Biomatrix Alpha.




➤ “Crushing stent” technique:

➤ Result after “Kissing balloon”:

TIMI 3 flow in both arteries with opening of the struts of the stent towards CX.

# Left main bifurcation treated with Biomatrix Alpha.



➤ “Crushing stent” technique:

➤ Optimization of the results:

“POTS” with two 3.0 NC balloon towards distal LM and prox LAD.






# Left main bifurcation treated with Biomatrix Alpha.

➤ “Crushing stent” technique:


➤ Optimization of the results:

IVUS to distal LM and prox LAD

Good aposition of the stent but underexpansion in LM.



# Left main bifurcation treated with Biomatrix Alpha.



➤ “Crushing stent”  
technique:

➤ Optimization of  
the results:

Final POTS with 4.0  
NC balloon to LM.



# Left main bifurcation treated with Biomatrix Alpha.

- Myocardial infarction represents a challenging scenario in PCI due to:
  - Clinical instability of patients.
  - &
  - Complex lesion: thrombus, invisible distal vessels..

# Left main bifurcation treated with Biomatrix Alpha.

- The presence of LM stenosis in a patient with STEMI always introduce a plus of complexity in the management of the case.
- The strategy in these cases is always related to clinical stability and operator skills.



# Left main bifurcation treated with Biomatrix Alpha.

- Percutaneous treatment of LM lesion has proven to be an effective strategy in elective patients (NOBLE trial, EXCEL Trial..).
- However, in the context of STEMI involving LM, the percutaneous treatment is probably the only chance.

# Left main bifurcation treated with Biomatrix Alpha.

- The present case is a challenging case due to the combination of STEMI plus LM involvement.
- *Biomatrix Alpha* has proven to be an excellent option in treating difficult lesion with complex techniques.



# Left main bifurcation treated with Biomatrix Alpha.

- Some of the remarkable characteristics of the stent demonstrated in this case are:
  - Excellent deployment capability in calcified and tortuous scenarios.
  - Strut structure that allows easy re-crossing after expansion.
    - Excellent overexpansion in big vessels.

# **Left main bifurcation treated with Biomatrix Alpha.**

- In conclusion, Biomatrix Alpha has demonstrated its adaptation to complex lesion, complex techniques, with excellent results.
- The drug and the polymer have demonstrated a safe profile in all scenarios where a DES is needed.